

Aspen Thoraco Lumbo Sacral Orthosis (Aspen TLSO Brace)

NHS Lothian /
South East Trauma Network



FRONT



BACK



Purpose of the brace

Following the injury to your spine the brace helps limit movement and provides stability while your injury heals.

Continue wearing the brace as instructed until this is reviewed, and the plan adjusted by your medical team.

You need to put the brace on lying flat:	Yes	No
You can put the brace on in sitting:	Yes	No
You are required to wear the brace at all times (Inc when sitting up):	Yes	No
You will need to sleep with the brace on:	Yes	No

You need to the wear the brace for a period of ____ weeks, until advised otherwise at your fracture clinic appointment.

Aspen Thoraco Lumbo Sacral Orthosis (Aspen TLSO Brace)

NHS Lothian /
South East Trauma Network



How to put the brace on:

1. Lying on your back, roll onto your side, aiming to keep your spine aligned; roll by holding your arms across your chest (hands to opposite shoulders) and bend your knees.



2. Position the back section so that it is centrally aligned. With your helper, ensure the bottom edge is in line with the base of your lower back.

3. They can press down on the mattress to feed the brace under your side.



4. Holding the brace in place roll onto your back. Ensure the brace is located centrally – the curve within the back piece should follow the natural curve of your lower spine.

5. Put the front section in place, this may already be attached to the left front strap, making sure it is centrally aligned with your waist and abdomen; then fold over the front straps, the right side will go on top/last.


6. At this point you may need to make small adjustments to ensure the brace is correctly placed.

Aspen Thoraco Lumbo Sacral Orthosis (Aspen TLSO Brace)


NHS Lothian /
South East Trauma Network



7. Now pull the cords directly up to tighten the brace until it feels snug, then fold them over the front and velcro them down.



8. Bring the shoulder straps over and click them into place.



9. The lower end of the brace should sit just above the groin crease and just above your buttocks at the back.

How to remove the brace

- 1.** Lie flat on your back, undo the Velcro cord tabs and then pull open the front section.
- 2.** Roll onto your side keeping the spine aligned (same technique as above).
- 3.** Remove the brace.
- 4.** Roll onto your back keeping the spine aligned.



General Advice

- 1.** The brace should be kept on at all times except when lying down - unless the consultant has advised you otherwise.
- 2.** The brace should be tight but comfortable.
- 3.** The brace should be worn over thin, preferably cotton clothing such as a T-Shirt or vest. You should avoid clothing with buttons or jewellery that would be pressed against your body by the brace.
- 4.** When sitting on a chair, the bottom of the brace at the back should be just sitting above the seat of the chair, if this is catching the chair, it will cause the brace to rise up.
- 5.** Skin checks should be carried out by a family member. Remove the brace for short periods during the day when lying down to allow your skin to breathe. If there is an area of broken skin, please get it reviewed by a practice nurse or GP and contact the Orthotics Department to have the brace reviewed.
- 6.** Use pain relief in the first few weeks to allow you to move comfortably. If high levels of pain persist, contact your GP.
- 7.** The brace can be cleaned by wiping it over with a damp cloth and leaving it to dry; do not place on a radiator or a fire to dry as this can damage the brace structure.



Daily Activities Advice

1. Bathing and showering are not permitted; your occupational therapist or physiotherapist can further discuss this with you.
2. Position changes and pacing yourself throughout the day can help manage your pain; it is best to do 'little and often'.
3. Avoid heavy lifting (no more than 5kg) and prolonged bending or twisting movements.
4. To put on socks and shoes, sit down and place your foot on your knee. You can consider a long handle shoehorn. When dressing your lower body, we recommend you sit down and bring your legs up bending at the hips and knees.
5. You can start light housework tasks as your pain allows. Kitchen tasks are best done at a counter top surface to avoid bending or overreaching.
6. Walking is recommended to aid your recovery. It will increase your activity levels and help reduce stiffness. Start with short distances and gradually increase as able.
7. Returning to work depends on your recovery and tasks involved in your job. Returning to heavy, manual tasks may take longer.
8. You are not able to drive as long as you are required to wear the brace.
9. Unless told otherwise follow up physiotherapy is initially not required. Please discuss the need for future physiotherapy at your fracture clinic appointment.

This leaflet is not meant to replace any discussion with your health care team, please let us know if you have any questions